



Cold Brew Black Tea																																			
Item Code: FPSPC024																																			
Uses: Tea, Drinks, Smoothie, Ice Cream, Dairy, Dessert		Benefits in Use: 100 % natural, cold Brew Extract																																	
Shelf Life Frozen: 1 Year	Brix: 8 – 12	Pack Volume: 50 mL																																	
Shelf Life Refrigerated: 21 days	pH: 5.10 - 5.15	Processing: Cold Pressed Cold Extraction																																	
		Ingredient Statement: Black Tea leaves	Regulatory: Allergen Free Non-GMO Made in Canada from Imported Ingredients																																
		Raw Material Country of Origin: China	GFSI Certification: FSSC 22000 Compliant																																
Microbiological Specifications:																																			
<table border="1"> <thead> <tr> <th>Description</th> <th>Acceptable Level</th> <th>Detection Limit</th> <th>Method/Test Used</th> </tr> </thead> <tbody> <tr> <td>ACC</td> <td><3000</td> <td><10 cfu/mL</td> <td>Petrifilm - MFHPB-18</td> </tr> <tr> <td>Yeast</td> <td><1000</td> <td><10 cfu/mL</td> <td>Petrifilm - MFHPB-22</td> </tr> <tr> <td>Mold</td> <td><300</td> <td><10 cfu/mL</td> <td>Petrifilm - MFHPB-22</td> </tr> <tr> <td><i>E. coli</i></td> <td><10</td> <td><10 cfu/mL</td> <td>Petrifilm - MFHPB-34</td> </tr> <tr> <td>Coliforms</td> <td><10</td> <td><10 cfu/mL</td> <td>Petrifilm - MFHPB-34</td> </tr> <tr> <td><i>Listeria sp</i></td> <td>Negative</td> <td>Negative</td> <td>MFLP-49</td> </tr> <tr> <td><i>Salmonella sp.</i></td> <td>Negative</td> <td>Negative</td> <td>MFLP-30</td> </tr> </tbody> </table>				Description	Acceptable Level	Detection Limit	Method/Test Used	ACC	<3000	<10 cfu/mL	Petrifilm - MFHPB-18	Yeast	<1000	<10 cfu/mL	Petrifilm - MFHPB-22	Mold	<300	<10 cfu/mL	Petrifilm - MFHPB-22	<i>E. coli</i>	<10	<10 cfu/mL	Petrifilm - MFHPB-34	Coliforms	<10	<10 cfu/mL	Petrifilm - MFHPB-34	<i>Listeria sp</i>	Negative	Negative	MFLP-49	<i>Salmonella sp.</i>	Negative	Negative	MFLP-30
Description	Acceptable Level	Detection Limit	Method/Test Used																																
ACC	<3000	<10 cfu/mL	Petrifilm - MFHPB-18																																
Yeast	<1000	<10 cfu/mL	Petrifilm - MFHPB-22																																
Mold	<300	<10 cfu/mL	Petrifilm - MFHPB-22																																
<i>E. coli</i>	<10	<10 cfu/mL	Petrifilm - MFHPB-34																																
Coliforms	<10	<10 cfu/mL	Petrifilm - MFHPB-34																																
<i>Listeria sp</i>	Negative	Negative	MFLP-49																																
<i>Salmonella sp.</i>	Negative	Negative	MFLP-30																																